



NEW TECHNOLOGY **INNOVATIVE** PSP SCANNER




Check the patient's image on the chair side quickly and easily



- ✔ Image preview on touch screen
- ✔ Image resolution of up to 17lp/mm provides image with stunning detail.
- ✔ Images are saved on built-in memory card and user can import images to PC when you want.
- ✔ Support fast response to customer complaint.
- ✔ Support multi-user practices.
- ✔ Support sleep mode for power saving.
- ✔ Light, compact and competitive price.
- ✔ Easy to install and use.




High quality **High resolution images**



- Flexible
- High patients comfort




🔍 Imaging plates




22 X 31mm



24 X 40mm



31 X 41mm



27 X 54mm

Make more accurate diagnosis with clear image quality with CRX-1000



CRX-1000 Specification

Supported IP	Size 0 / 1 / 2 / 3
Resolution(lp/mm)	17
Display	4" Touch Screen
Interface	Ethernet
Multi-user Practices	Yes
Sleep Mode	Yes
Weight(kg)	3.3
Dimensions(H x W x Dmm)	233 x 141 x 191.5